



***Nanosecond Fluorescence Lifetime
Imaging using a directly gated interline
transfer CCD detector
by
Photonic Research Systems Ltd***


- **Sensor: 2/3" format monochrome interline CCD sensor, peltier cooled to -18°C**
- **Light Source: Diode pumped, frequency tripled solid-state pulsed laser with active Q-Switching. Pulse repetition rate was 5kHz, wavelength 355nm, pulse width at half maximum was 1.2 nsec, laser jitter was <1nsec.**
- **Samples: Two solutions of fluorescein were prepared one having twice the emission intensity of the other. The brighter sample was then quenched with potassium iodide until its intensity matched that of the lower intensity sample**


Cuvettes of fluorescein illuminated by pulsed UV laser



One cuvette holds fluorescein in buffer (4ns) while the other has a sample of twice the concentration but quenched with potassium iodide to match the dilute sample (2ns). This image was taken with the *Imagex* camera in steady state mode under laser illumination and with some ambient light.

Steady-State Fluorescence image of Fluorescein samples with no background light





2ns delay


4ns delay


6ns delay

Delay following laser pulse


The next slide shows the resultant decay-weighted image when the appropriate pair of time-resolved images are divided. A psuedo-colour look-up table is used to show the lifetime differences.

Note the very high signal-to-noise ratio, which is evident even in the uniform colour of the pixels representing very weak signals from reflections at the meniscus and base of the cuvettes





The Imagex Matrix Dialog allows the user to select 2 source images and adjust their contrast settings. The images are then combined into a single 'Matrix' image



***Combining the Intensity and Flim Images
in this way helps to suppress artifacts created by performing
ratio operations on areas of very low intensity.***



quenched unquenched

The image shows two horizontal fluorescence bands. The left band is labeled 'quenched' and appears as a faint purple line. The right band is labeled 'unquenched' and appears as a bright yellow-green line. The background is dark, and the bands are centered horizontally.

- ***Investigation of smaller format and higher resolution CCD sensors***
- ***Investigate Improved DPSS Light Source:
0.8 nsec FWHM with jitter <0.5nsec now available.***
- ***Mount System on Microscope for real biological applications!***